

# Human GLP1R Full Length Protein (L260F), Flag,His Tag (Detergent)

Catalog # GLR-H52D4



BIOSYSTEMS  
**Acro**

## Synonym

GLP-1-R, GLP-1R, glucagon like peptide 1 receptor

## Source

Human GLP1R Full Length Protein (L260F), Flag,His Tag(GLR-H52D4) is expressed from human 293 cells (HEK293). It contains AA Arg 24 - Ser 432 (Accession # [P43220-1](#)(L260F)).

Predicted N-terminus: Asp

## Molecular Characterization



This protein carries flag tag at the N-terminus and polyhistidine tag at the C-terminus.

The protein has a calculated MW of 66.5 kDa.

## Endotoxin

Less than 1.0 EU per  $\mu\text{g}$  by the LAL method.

\*The DDM/CHS buffer (Cat. No. [DC-11](#)) is sold separately and not included in protein, you can follow [this link](#) for product information.

## Formulation

*This product is not suitable for cell based experiments due to cytotoxicity of DDM.*

*DDM and CHS are INDISPENSABLE to keep membrane protein soluble and active, under no circumstance should you remove DDM and CHS.*

*DDM/CHS buffer (DC-11) is sold separately and not included in protein, and please contact us if you need the buffer.*

*If glycerol is not compatible to your application, remove glycerol just before immediate experiment, and NEVER store glycerol-free protein solution.*

Supplied as 0.2  $\mu\text{m}$  filtered solution in 50 mM HEPES, 150 mM NaCl, DDM, CHS, pH7.5 with glycerol as protectant.

Contact us for customized product form or formulation.

## Shipping

*This product is supplied and shipped with dry ice, please inquire the shipping cost.*


## Storage

*Please avoid repeated freeze-thaw cycles.*

This product is stable after storage at:

- The product MUST be stored at  $-70^{\circ}\text{C}$  or lower upon receipt;
- $-70^{\circ}\text{C}$  for 3 months under sterile conditions.

## Bioactivity-ELISA



Immobilized Human GLP1R Full Length Protein, Flag,His Tag (Cat. No. GLR-H52D4) at 5  $\mu\text{g}/\text{mL}$  (100  $\mu\text{L}/\text{well}$ ) on a Nickel Coated plate (Cat. No. SP-16) can bind GLP1R agonist with a linear range of 0.5-16 ng/mL (QC tested).


## Bioactivity-SPR

Discounts, Gifts,  
and more!



[www.acrobiosystems.com](http://www.acrobiosystems.com)

8/7/2024



GLP1R agonist captured on Protein A Chip can bind Human GLP1R Full Length Protein, Flag, His Tag (Cat. No. GLR-H52D4) with an affinity constant of 2.74  $\mu$ M as determined in a SPR assay (in presence of DDM and CHS) (Biacore 8K) (Routinely tested).

## Background

This gene encodes a 7-transmembrane protein that functions as a receptor for glucagon-like peptide 1 (GLP-1) hormone, which stimulates glucose-induced insulin secretion. This receptor, which functions at the cell surface, becomes internalized in response to GLP-1 and GLP-1 analogs, and it plays an important role in the signaling cascades leading to insulin secretion. It also displays neuroprotective effects in animal models. Polymorphisms in this gene are associated with diabetes. The protein is an important drug target for the treatment of type 2 diabetes and stroke.

## Clinical and Translational Updates

Discounts, Gifts,  
and more!

