



Synonym

ILDR2,C1orf32

Source

Human ILDR2, Fc Tag(IL2-H5259) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Glu 186 (Accession # [Q71H61-1](#)).

Predicted N-terminus: Leu 21

Molecular Characterization

| | |
|-------------------------------------|---------------------------------|
| ILDR2(Leu 21 - Glu 186) Q71H61-1 | Fc(Pro 100 - Lys 330) P01857 |
|-------------------------------------|---------------------------------|

This protein carries a human IgG1 Fc tag at the C-terminus.

The protein has a calculated MW of 45.1 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>95% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Human ILDR2, Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

Background

ILDR2 belongs to the immunoglobulin superfamily. May be involved in lipid homeostasis and ER stress pathways. ILDR2 as a novel B7-like protein with robust T cell inhibitory activity, expressed in immune cells and in immune-privileged and inflamed tissues. A fusion protein, consisting of ILDR2 extracellular domain with an Fc fragment, that binds to a putative counterpart on activated T cells showed a beneficial effect in the collagen-induced arthritis model and abrogated the production of proinflammatory cytokines and chemokines in autologous synovial-like cocultures of macrophages and cytokine-stimulated T cells. Collectively, these

Discounts, Gifts,
and more!



Human ILDR2 Protein, Fc Tag

Catalog # IL2-H5259



findings point to ILDR2 as a novel negative regulator for T cells, with potential roles in the development of immune-related diseases, including autoimmunity and cancer.

Clinical and Translational Updates

Discounts, Gifts,
and more!



www.acrobiosystems.com

6/20/2024